

Product datasheet for **MR206877L3V**

Sept4 (BC055101) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Sept4 (BC055101) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Sept4 |
| Synonyms: | ARTS; Bh5; OTTMUSG00000001265; Pnutl2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | BC055101 |
| ORF Size: | 1293 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR206877). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | BC055101.1 |
| RefSeq Size: | 1680 bp |
| RefSeq ORF: | 1295 bp |
| Locus ID: | 18952 |
| Cytogenetics: | 11 C |



[View online »](#)

Gene Summary:

Filament-forming cytoskeletal GTPase. Forms a filamentous structure with SEPTIN12, SEPTIN6, SEPTIN2 and probably SEPTIN4 at the sperm annulus which is required for the structural integrity and motility of the sperm tail during postmeiotic differentiation (By similarity). May play a role in cytokinesis (Potential). May play a role in platelet secretion (By similarity).[UniProtKB/Swiss-Prot Function]