

Product datasheet for **MR206897L3V**

Fgl2 (NM_008013) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgl2 (NM_008013) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgl2
Synonyms:	AI385601; musfiblp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008013
ORF Size:	1299 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206897).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008013.2
RefSeq Size:	3769 bp
RefSeq ORF:	1299 bp
Locus ID:	14190
UniProt ID:	P12804
Cytogenetics:	5 9.83 cM
Gene Summary:	Converts prothrombin to thrombin.[UniProtKB/Swiss-Prot Function]


[View online »](#)