

## Product datasheet for **MR206947L3V**

### Stambpl1 (NM\_029682) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Stambpl1 (NM_029682) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stambpl1
Synonyms:	1700095N21Rik; 8230401J17Rik; ALMalpha; AMSH-FP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029682
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206947).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029682.4</a>
RefSeq Size:	1978 bp
RefSeq ORF:	1311 bp
Locus ID:	76630
UniProt ID:	<a href="#">Q76N33</a>
Cytogenetics:	19 C1
Gene Summary:	Zinc metalloprotease that specifically cleaves 'Lys-63'-linked polyubiquitin chains. Does not cleave 'Lys-48'-linked polyubiquitin chains (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)