

## Product datasheet for **MR206950L3V**

### Fgg (NM\_133862) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgg (NM_133862) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgg
Synonyms:	3010002H13Rik; AI256424
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133862
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206950).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133862.1</a>
RefSeq Size:	1536 bp
RefSeq ORF:	1311 bp
Locus ID:	99571
UniProt ID:	<a href="#">Q8VCM7</a>
Cytogenetics:	3 36.94 cM



[View online »](#)

**Gene Summary:**

This gene encodes the gamma subunit of the coagulation factor fibrinogen, which is a component of the blood clot. Mice lacking the encoded protein did not possess detectable amounts of plasma fibrinogen. Pregnant mice lacking the encoded protein die due to heavy bleeding during delivery. This gene is located adjacent to the genes encoding fibrinogen alpha and gamma subunits on chromosome 3. Alternate splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Nov 2015]