


Product datasheet for MR206953L3

Stk11 (NM_011492) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Stk11 (NM_011492) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Stk11
Synonyms:	AA408040; Lkb1; Par-4; R75140
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206953).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.

ACCN:	NM_011492
ORF Size:	1308 bp



OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_011492.3
RefSeq Size:	1643 bp
RefSeq ORF:	1311 bp
Locus ID:	20869
UniProt ID:	Q9WTK7
Cytogenetics:	10 C1
Gene Summary:	This gene encodes a member of the serine/threonine kinase family. The encoded protein, a known tumor suppressor, activates (via phosphorylation) adenine monophosphate-activated protein kinase (AMPK) and AMPK-related kinase proteins. This upstream regulation of the AMPK pathway is thought to regulate a number of different processes, including cell metabolism, cell polarity, apoptosis and DNA damage response. Mutations in a similar gene in human have been associated with Peutz-Jeghers syndrome. Alternative splicing results in multiple transcript variants, including the S isoform which plays a potential role in spermiogenesis. [provided by RefSeq, Sep 2014]

Product images:



Circular map for MR206953L3