

## Product datasheet for **MR206972L3V**

### Plin3 (NM\_025836) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Plin3 (NM_025836) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plin3
Synonyms:	1300012C15Rik; M6prbp1; Tip47
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025836
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206972).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025836.3</a> , <a href="#">NP_080112.1</a>
RefSeq Size:	2159 bp
RefSeq ORF:	1314 bp
Locus ID:	66905
UniProt ID:	<a href="#">Q9DBG5</a>
Cytogenetics:	17 D
Gene Summary:	Required for the transport of mannose 6-phosphate receptors (MPR) from endosomes to the trans-Golgi network.[UniProtKB/Swiss-Prot Function]



[View online »](#)