

Product datasheet for MR206977L4

Alg1 (BC011281) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Alg1 (BC011281) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Alg1
Synonyms:	HMT1, HMAT1, MGC18946
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206977).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	BC011281
ORF Size:	1316 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC011281
RefSeq Size:	1677 bp
RefSeq ORF:	1316 bp
Locus ID:	208211
Cytogenetics:	16 A1
Gene Summary:	Participates in the formation of the lipid-linked precursor oligosaccharide for N-glycosylation. Involved in assembling the dolichol-pyrophosphate-GlcNAc(2)-Man(5) intermediate on the cytoplasmic surface of the ER (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR206977L4