

## Product datasheet for **MR206987L3V**

### Alg3 (NM\_145939) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Alg3 (NM_145939) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Alg3
Synonyms:	D16Ertd36e; not
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145939
ORF Size:	1317 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206987).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145939.1</a> , <a href="#">NP_666051.1</a>
RefSeq Size:	1425 bp
RefSeq ORF:	1317 bp
Locus ID:	208624
UniProt ID:	<a href="#">Q8K2A8</a>
Cytogenetics:	16 12.48 cM
Gene Summary:	Adds the first Dol-P-Man derived mannose in an alpha-1,3 linkage to Man5GlcNAc2-PP-Dol. [UniProtKB/Swiss-Prot Function]



[View online »](#)