

## Product datasheet for **MR206991L4V**

### Smo (BC048091) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Smo (BC048091) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Smo
Synonyms:	D13Mgi8, Smoh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC048091
ORF Size:	1317 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206991).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC048091</a> , <a href="#">AAH48091</a>
RefSeq Size:	1747 bp
RefSeq ORF:	1319 bp
Locus ID:	319757
Cytogenetics:	6 12.36 cM



[View online »](#)

**Gene Summary:**

G protein-coupled receptor that probably associates with the patched protein (PTCH) to transduce the hedgehog's proteins signal. Binding of sonic hedgehog (SHH) to its receptor patched is thought to prevent normal inhibition by patched of smoothened (SMO) (By similarity). Required for the accumulation of KIF7, GLI2 and GLI3 in the cilia. Interacts with DLG5 at the ciliary base to induce the accumulation of KIF7 and GLI2 at the ciliary tip for GLI2 activation (PubMed:25644602).[UniProtKB/Swiss-Prot Function]