

Product datasheet for **MR207010**

Cpne8 (BC017626) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cpne8 (BC017626) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cpne8
Synonyms:	1200003E11Rik; 1500031E20Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR207010 representing BC017626
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGTTGATTCTAAGAGTCCCACTTGTCCTCAACATAGGAATTCAGGCAGGAAATGTGGACAATCATAC
 TGACGGCAGAGGAGCTAACTGCTGCAGGGATGCTGTCTTGATGCAGTTCTGTGCAACAAATTGGATAA
 GAAGGATTTCTTTGGGAAGTCAGATCCTTTCTGGTATTTTATCGAAGCAATGAGGATGGAAGCTTTACA
 ATTTGCCACAAGACAGAAGTGGTCAAAACACTCTAAATCCTGTGTGGCAAGCGTTCAAGATCTCAGTGC
 GAGCATTGTGCAATGGGGACTATGACAGAACAATCAAGGTAGAAGGTACGACTGGGACCGAGACGGGAG
 TCATGATTTTCATTGGAGAATTCACAACAAGCTATCGGGAGCTCGCAAGAGGGCAGTCACAATTCAACGTC
 TATGAGGTGGTGAATCCCAAAAAGAAAGGAAAAAAGAAAAATACACCAATTCAGGAACAGTAACCTTAT
 TGTCTTCTTGGTGGAAACAGAAGTCTCCTTCTGGACTATATTAAGGAGGAACGCAGATCAACTTCAC
 GGTGGCCATTGATTTTACAGCATCAAATGGAAACCTGCCAGCCAATCTCTCCACTATATGAACCTT
 TACCAACTGAATGCCTATGGAATGGCACTGAAGGCTGTAGGAGAAATCGTTCAAGATTATGACAGTGATA
 AGATGTTCCCGCCCTGGGTTTTGGTGCAGGCTGCCTCCAGATGGAAGGATATCTCACGAGTTGCCCT
 GAATGGCAACCTCAGAATCCTTACTGTGATGGCATTGAAGGTGTGATGGAAGCTTACTATAGAAGTCTG
 AAGTCTGTCCAGTTGTATGGGCCAACCACTTTGCTCCCGTAATAAACACGTAGCAAGATACGCATCCT
 CTGTAAAGATGGCTCCAGTACTTTGTCTTCTCATTGTGACGGACGGCGTCATCTCAGATATGGCTCA
 AACTAAAGATCCATAGTTAATGCCTCAAAGCTTCCCATGTCAATAATAAGTAGGTGTGGACCTGCA
 GAATTTGATGCAATGGTTGAATTGGATGGAGATGATGTAAGAGTTTCTCCAGAGGAAAATATGCTGAAA
 GAGATATTGTACAGTTTGTGCCGTTTCAGGATTATATTGACAGAAGTGGAAACCATATCAGCATGGC
 GAGACTGGCCAAAGATGTCTTAGCTGAGATCCAGAGCAGTTTCTTTCCTACATGAGAGCCAGAGGAATC
 AAGCCATACCGGCGCCCCCGCGGTACACCCCTCTACACACGTGCTGCAGACGCAATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR207010 representing BC017626
 Red=Cloning site Green=Tags(s)

MLILRVPTCPNIGIPGRKCGTIILTAEELNCCRDVLMQFCANKLDKKDFFGKSDPFLVFYRSNEDGSFT
 ICHKTEVVKNLNPVWQAFKISVRALCNGDYDRTIKVEVYDWRDGSDFIGFTTSYRELARGQSQFNV
 YEYVNPKKKGGKKKYTNSTVTLSSFLVETEVSFLDYIKGGTQINFVAIDFTASNGNPAQPTSLHYMNP
 YQLNAYGMALKAVGEIVQDYSDKMFALGFGAKLPDGRISHEFALNGNPQNPYCDGIEGVMEAYYRSL
 KSVQLYGPTNFAPVINHVARYASSVKDGSQYFVLLIVTDGVIDMAQTKEIVNASKLPMSIIIVGVGPA
 EFDAMVELDGDVVRSSRGKYAERDQVFPFRDYIDRSNHLISMARLAKDVLAEIPEQFLSYMRARGI
 KPSPAPPYPPTPHVLQTQI

TRTRPLE**QKLISEEDLA**AND**ILDYKDDDDKV**

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: BC017626

ORF Size: 1320 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC017626.1](#)


RefSeq Size: 3231 bp

RefSeq ORF: 1322 bp

Locus ID: 66871

Cytogenetics:	15 E3
MW:	118.5 kDa
Gene Summary:	This gene encodes a member of the copine family of highly conserved, calcium-dependent phospholipid binding proteins. The encoded protein has two characteristic C2 domains and a VWFA domain and may play a role in membrane trafficking. A related pseudogene is found on chromosome 8. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR207010