

Product datasheet for **MR207034L3V**

Gpr177 (BC018381) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpr177 (BC018381) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gpr177
Synonyms:	5031439A09Rik; AI173978; AI987742; Gpr177
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC018381
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207034).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC018381 , AAH18381
RefSeq Size:	2188 bp
RefSeq ORF:	1325 bp
Locus ID:	68151
Cytogenetics:	3 H4



[View online »](#)

Gene Summary:

Regulates Wnt proteins sorting and secretion in a feedback regulatory mechanism. This reciprocal interaction plays a key role in the regulation of expression, subcellular location, binding and organelle-specific association of Wnt proteins. Plays also an important role in establishment of the anterior-posterior body axis formation during development.
[UniProtKB/Swiss-Prot Function]