

Product datasheet for MR207060

Ehd2 (BC027084) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ehd2 (BC027084) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ehd2
Synonyms: MGC25606, MGC38650
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR207060 representing BC027084
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCATGGCGAGACAGAGGGCACCGTGCCCGGAAATGCTCTCGTTGTGGACCCAGAGAAGCCCTCCGCA
 AACTCAACCCCTTTGGAAACACTTTTCTCAACAGGTTTATGTGCGCCCAACTCCCAACAGGTCCTGGA
 GAGCATCAGCATCATCGACACGCTGGCATCCTATCAGGTGCAAAACAGAGAGTGAGCCGTGGCTACGAC
 TTCCCAGCTGTCTGCGCTGTTTGTGCTGAACGCGTGGACTCATCATCTTGCTCTTTGATGCCACAAGC
 TGGAGATTTCCGATGAATTCTCGGAGGCCATCGGTGCCTTGCCTGGCCACGAGGATAAGATCCGCGTGGT
 ACTCAACAAGGCAGACATGGTGGAGACACAGCAGCTGATGCGCGTCTACGGGGCGCTCATGTGGGCGCTG
 GGCAAGGTGGTGGGCACGCTGAGGTCCTTCGCGTCTACATTGGCTCCTTCTGGTCCCAGCCGCTCCTTG
 TGCCGGACAACCGGCTCTTTTCGAGCTGGAGGAACAGGACCTTCCGTGATATCCAAGGCCTACCGCG
 CCATGCGGCGCTGCGCAAGCTCAACGACCTAGTGAAACGGGCGCGGCTGGTGCAGGACACGCTTACATC
 ATCAGCTACCTGAAGAAAGAGATGCCACGCTGTTTCGGGAAGGAAAACAAGAAGAAGCAGCTGATCCTGA
 AGTGCCTGTCACTTTTGGGAAGATCCAGCTGGAGCACCATCTCCCCTGGAGACTTTCTGACTGCCA
 GAAGATGCAGGAGCTGTTGATGGCACATGACTTCACCAAGTTCCTACTCGCTAAAGCCGAAGCTGCTGGAG
 GCGCTAGATGACATGCTGGCCCAAGACATTGCCAAGCTTATGCCGCTGCTTCGGCAGGAGGAGCTGGAGA
 GTGTGAGGCTGGAGTGCAGGGCGGTGCGTTTGGAGGACCCGCGATGGGCCCCTTCGTGGAGCGGGGCC
 TGATGAAGCAATAGAGGATGGTGGAGGGGCTCAGAAGATGATGCCGAATGGTGGTACCAAGGACAAG
 TCCAAATACGATGAAATCTTCTACAACCTAGCTCCTGCTGATGGCAAGCTGAGTGGCTCTAAGGCCAAA
 CCTGGATGGTGGGCACAAAGCTACCAACTCTGTGCTGGGAGAACTGGAAGCTGAGCGATGTCGACCG
 AGACGGCATGCTGGATGATGAGGAGTTTGCCTTGCCAGCCACCTCATTGAGGCTAAGCTCGAAGGCCAT
 GGTCTACCCACCAACCTACCCCGCCGCTGGTCCGCTTCTAAGCGACGACAGAAGGGCTCTGTGAA

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTAA




Protein Sequence: >MR207060 representing BC027084
Red=Cloning site Green=Tags(s)

MHGETEGTVPGNALVVDPEKPFRKLNPFNGNTFLNRFMCAQLPNQVLESISIIDTPGILSGAKQRVSRGYD
 FPAVLRWFAERVDLIILLFDAHKLEISDEFSEATGALRGHEDKIRVVLNKADMVETQQLMRVYGALMWAL
 GKVVGTPEVLRVYIGSFWSQPLLVPDNRRLELEEQDLFRDIQGLPRHAALRKLNDLVKRARLVRVHAYI
 ISYLKEMPTVFGKENKKKQLILKLPVIFAKIQLEHHISPGDFPDCQKMQELLMAHDFTKFHSLKPKLLE
 ALDDMLAQDIAKLMPLLRQEELESVEAGVQGGAFEGTRMGPFVERGPDEAIEDGEEGSEDDAEWVTKDK
 SKYDEIFYNLAPADGKLSGSKAKTWMVGTKLPSVLRGRIWKLSDVDRDGLMDDEEFALASHLIEAKLEGH
 GLPTNLPRLVPPSKRRQKGSAAE


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: BC027084

ORF Size:	1329 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC027084.1
RefSeq Size:	2717 bp
RefSeq ORF:	1331 bp
Locus ID:	259300
Cytogenetics:	7 8.65 cM
MW:	50 kDa
Gene Summary:	ATP- and membrane-binding protein that controls membrane reorganization/tubulation upon ATP hydrolysis (PubMed:24508342). Plays a role in membrane trafficking between the plasma membrane and endosomes. Important for the internalization of GLUT4 (PubMed:14676205). Required for fusion of myoblasts to skeletal muscle myotubes. Required for normal translocation of FER1L5 to the plasma membrane (PubMed:18502764, PubMed:21177873). Regulates the equilibrium between cell surface-associated and cell surface-dissociated caveolae by constraining caveolae at the cell membrane (By similarity). [UniProtKB/Swiss-Prot Function]