

## Product datasheet for **MR207074L3V**

### Efemp2 (NM\_001164352) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Efemp2 (NM_001164352) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Efemp2
Synonyms:	0610011K11Rik; Fbln4; MBP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164352
ORF Size:	1332 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207074).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164352.1</a> , <a href="#">NP_001157824.1</a>
RefSeq Size:	1802 bp
RefSeq ORF:	1389 bp
Locus ID:	58859
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

Plays a crucial role in elastic fiber formation in tissue, and in the formation of ultrastructural connections between elastic laminae and smooth muscle cells in the aorta, therefore participates in terminal differentiation and maturation of smooth muscle cell (SMC) and in the mechanical properties and wall integrity maintenance of the aorta (PubMed:16478991, PubMed:19855011, PubMed:20019329, PubMed:26486174, PubMed:26711913, PubMed:28508064). In addition, is involved in the control of collagen fibril assembly in tissue through proteolytic activation of LOX leading to cross-linking of collagen and elastin (PubMed:26690653, PubMed:26711913, PubMed:26220971, PubMed:26178373). Also promotes ELN coacervation and participates in the deposition of ELN coacervates on to microfibrils but also regulates ELN cross-linking through LOX interaction (PubMed:17324935). Moreover adheres to the cells through heparin binding in a calcium-dependent manner and regulates vascular smooth muscle cells proliferation through angiotensin signaling (PubMed:23636094).[UniProtKB/Swiss-Prot Function]