

Product datasheet for **MR207130L4V**

Bag5 (NM_027404) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bag5 (NM_027404) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bag5
Synonyms:	1600025G07Rik; 4930405J06Rik; Bag-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027404
ORF Size:	1344 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207130).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027404.1
RefSeq Size:	2150 bp
RefSeq ORF:	1344 bp
Locus ID:	70369
UniProt ID:	Q8CI32
Cytogenetics:	12 F1
Gene Summary:	May function as a nucleotide exchange factor for HSP/HSP70, promoting ADP release, and activating Hsp70-mediated refolding. Inhibits both auto-ubiquitination of PRKN and ubiquitination of target proteins by PRKN (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)