

Product datasheet for **MR207142L4V**

Ispd (NM_178629) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ispd (NM_178629) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ispd
Synonyms:	4930579E17Rik; AV040780
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178629
ORF Size:	1344 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207142).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178629.5 , NP_848744.2
RefSeq Size:	2818 bp
RefSeq ORF:	1344 bp
Locus ID:	75847
UniProt ID:	Q5RJG7
Cytogenetics:	12 A3



[View online »](#)

Gene Summary:

Cytidylyltransferase required for protein O-linked mannosylation (By similarity). Catalyzes the formation of CDP-ribitol nucleotide sugar from D-ribitol 5-phosphate (By similarity). CDP-ribitol is a substrate of FKTN during the biosynthesis of the phosphorylated O-mannosyl trisaccharide (N-acetylgalactosamine-beta-3-N-acetylglucosamine-beta-4-(phosphate-6-)mannose), a carbohydrate structure present in alpha-dystroglycan (DAG1), which is required for binding laminin G-like domain-containing extracellular proteins with high affinity (By similarity). Shows activity toward other pentose phosphate sugars and mediates formation of CDP-ribulose or CDP-ribose using CTP and ribulose-5-phosphate or ribose-5-phosphate, respectively (By similarity). Not Involved in dolichol production (By similarity). [UniProtKB/Swiss-Prot Function]