

Product datasheet for **MR207173L4V**

Tuba8 (NM_017379) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tuba8 (NM_017379) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tuba8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017379
ORF Size:	1350 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017379.1
RefSeq Size:	1432 bp
RefSeq ORF:	1350 bp
Locus ID:	53857
UniProt ID:	Q9JIZ2
Cytogenetics:	6 F1
Gene Summary:	This gene encodes a member of the tubulin family of cytoskeletal proteins that form the integral component of microtubules. This gene is preferentially expressed in heart, skeletal muscle and testis. The encoded protein may play a role in liver tumorigenesis in mice. [provided by RefSeq, Sep 2015]



[View online »](#)