

Product datasheet for **MR207188L4V**

Lancl2 (BC016072) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lancl2 (BC016072) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lancl2
Synonyms:	1700003F10Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC016072
ORF Size:	1350 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207188).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC016072 , AAH16072
RefSeq Size:	2204 bp
RefSeq ORF:	1352 bp
Locus ID:	71835
Cytogenetics:	6 B3
Gene Summary:	Necessary for abscisic acid (ABA) binding on the cell membrane and activation of the ABA signaling pathway in granulocytes.[UniProtKB/Swiss-Prot Function]



[View online »](#)