

Product datasheet for **MR207224L3V**

Slc7a4 (BC016100) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc7a4 (BC016100) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc7a4
Synonyms:	MGC27672
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC016100
ORF Size:	1356 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207224).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC016100 , AAH16100
RefSeq Size:	1921 bp
RefSeq ORF:	1358 bp
Locus ID:	224022
Cytogenetics:	16 A3
Gene Summary:	Involved in the transport of the cationic amino acids (arginine, lysine and ornithine). [UniProtKB/Swiss-Prot Function]



[View online »](#)