

Product datasheet for **MR207233L4V**

Vrk3 (NM_133945) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vrk3 (NM_133945) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Vrk3
Synonyms:	AI428238
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133945
ORF Size:	1362 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207233).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133945.1 , NP_598706.1
RefSeq Size:	1790 bp
RefSeq ORF:	1362 bp
Locus ID:	101568
UniProt ID:	Q8K3G5
Cytogenetics:	7 B3
Gene Summary:	Inactive kinase that suppresses ERK activity by promoting phosphatase activity of DUSP3 which specifically dephosphorylates and inactivates ERK in the nucleus.[UniProtKB/Swiss-Prot Function]



[View online »](#)