

Product datasheet for **MR207244L3V**

Pfkfb4 (BC057594) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pfkfb4 (BC057594) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pfkfb4
Synonyms:	C230090D14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC057594
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207244).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC057594 , AAH57594
RefSeq Size:	3329 bp
RefSeq ORF:	1361 bp
Locus ID:	270198
Cytogenetics:	9 F2
Gene Summary:	Synthesis and degradation of fructose 2,6-bisphosphate.[UniProtKB/Swiss-Prot Function]



[View online »](#)