

## Product datasheet for **MR207247L3V**

### Casp9 (NM\_015733) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Casp9 (NM_015733) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Casp9
Synonyms:	A1115399; APAF-3; AW493809; Casp; CASP-9; Caspase-9; ICE-; ICE-LAP6; Mch6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015733
ORF Size:	1365 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015733.4</a>
RefSeq Size:	3899 bp
RefSeq ORF:	1365 bp
Locus ID:	12371
UniProt ID:	<a href="#">Q8C3Q9</a>
Cytogenetics:	4 D3



[View online »](#)

**Gene Summary:**

This gene is part of a family of caspases, aspartate-specific cysteine proteases well studied for their involvement in immune and apoptosis signaling. This protein, the initiator caspase, is activated after cytochrome c release from mitochondria and targets downstream effectors. In mouse, deficiency of this gene can cause perinatal lethality. This protein may have a role in normal brain development. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Apr 2013]