

Product datasheet for **MR207286L4V**

Smpdl3b (NM_133888) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Smpdl3b (NM_133888) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Smpdl3b
Synonyms:	1110054A24Rik; Asml3b; AU045240
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133888
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207286).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133888.2
RefSeq Size:	1925 bp
RefSeq ORF:	1371 bp
Locus ID:	100340
UniProt ID:	P58242
Cytogenetics:	4 D2.3



[View online »](#)

Gene Summary:

Lipid-modulating phosphodiesterase. Active on the surface of macrophages and dendritic cells and strongly influences macrophage lipid composition and membrane fluidity (PubMed:26095358). Acts as a negative regulator of Toll-like receptor signaling (PubMed:26095358, PubMed:27687724). Has in vitro phosphodiesterase activity, but the physiological substrate is unknown (PubMed:26095358, PubMed:27687724). Lacks activity with phosphocholine-containing lipids, but can cleave CDP-choline, and can release phosphate from ATP and ADP (in vitro) (PubMed:27687724).[UniProtKB/Swiss-Prot Function]