

Product datasheet for **MR207287L3V**

Rnh1 (NM_145135) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnh1 (NM_145135) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnh1
Synonyms:	AW546468; C80305; RNH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145135
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207287).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145135.2 , NP_660117.2
RefSeq Size:	1807 bp
RefSeq ORF:	1371 bp
Locus ID:	107702
UniProt ID:	Q91VI7
Cytogenetics:	7 F5
Gene Summary:	Ribonuclease inhibitor which inhibits RNASE1, RNASE2 and ANG. May play a role in redox homeostasis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)