

## Product datasheet for **MR207289L4V**

### Enpp4 (NM\_199016) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Enpp4 (NM_199016) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Enpp4
Synonyms:	4933413N07Rik; AA986363; AI195364; E-NPP 4; Gm90; mKIAA0879
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199016
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207289).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199016.2</a> , <a href="#">NP_950181.2</a>
RefSeq Size:	4558 bp
RefSeq ORF:	1371 bp
Locus ID:	224794
UniProt ID:	<a href="#">Q8BTJ4</a>
Cytogenetics:	17 B3



[View online »](#)

**Gene Summary:**

Hydrolyzes extracellular Ap3A into AMP and ADP, and Ap4A into AMP and ATP. Ap3A and Ap4A are diadenosine polyphosphates thought to induce proliferation of vascular smooth muscle cells. Acts as a procoagulant, mediating platelet aggregation at the site of nascent thrombus via release of ADP from Ap3A and activation of ADP receptors (By similarity). [UniProtKB/Swiss-Prot Function]