

Product datasheet for **MR207299L3V**

Angptl6 (NM_145154) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Angptl6 (NM_145154) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Angptl6
Synonyms:	6330404E11Rik; AGF; ARP3; Arp5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145154
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207299).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145154.2
RefSeq Size:	1564 bp
RefSeq ORF:	1374 bp
Locus ID:	70726
UniProt ID:	Q8R0Z6
Cytogenetics:	9 A3



[View online »](#)

Gene Summary:

May play a role in the wound healing process. May promote epidermal proliferation, remodeling and regeneration. May promote the chemotactic activity of endothelial cells and induce neovascularization. May counteract high-fat diet-induced obesity and related insulin resistance through increased energy expenditure.[UniProtKB/Swiss-Prot Function]