

Product datasheet for **MR207330L4V**

Snx8 (NM_172277) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snx8 (NM_172277) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snx8
Synonyms:	B130023O14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172277
ORF Size:	1380 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172277.1
RefSeq Size:	2652 bp
RefSeq ORF:	1380 bp
Locus ID:	231834
UniProt ID:	Q8CFD4
Cytogenetics:	5 G2
Gene Summary:	May be involved in several stages of intracellular trafficking. May play a role in intracellular protein transport from early endosomes to the trans-Golgi network (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)