

Product datasheet for **MR207333L3V**

Ggnbp1 (BC054732) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ggnbp1 (BC054732) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ggnbp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC054732
ORF Size:	1377 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207333).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC054732
RefSeq Size:	5028 bp
RefSeq ORF:	1379 bp
Locus ID:	70772
Cytogenetics:	17 A3.3
Gene Summary:	Induces mitochondrial fragmentation, possibly by promoting DNM1L-dependent fission and may play a role in mitochondrial morphogenesis during spermatogenesis.[UniProtKB/Swiss-Prot Function]


[View online »](#)