

## Product datasheet for **MR207339L4V**

### Chrm1 (NM\_001112697) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Chrm1 (NM_001112697) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Chrm1
Synonyms:	Chrm-1; M1; M1R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001112697
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207339).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001112697.1</a> , <a href="#">NP_001106167.1</a>
RefSeq Size:	6301 bp
RefSeq ORF:	1383 bp
Locus ID:	12669
UniProt ID:	<a href="#">P12657</a>
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

The muscarinic acetylcholine receptor mediates various cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels through the action of G proteins. Primary transducing effect is Pi turnover.[UniProtKB/Swiss-Prot Function]