

Product datasheet for **MR207349L3V**

Slc30a6 (NM_144798) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc30a6 (NM_144798) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc30a6
Synonyms:	9530029F08Rik; ZnT6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144798
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207349).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144798.2 , NP_659047.2
RefSeq Size:	2156 bp
RefSeq ORF:	1383 bp
Locus ID:	210148
UniProt ID:	Q8BJM5
Cytogenetics:	17 E2
Gene Summary:	Zinc-efflux transporter which allocates the cytoplasmic zinc to the trans-Golgi network (TGN) as well as the vesicular compartment.[UniProtKB/Swiss-Prot Function]



[View online »](#)