

Product datasheet for **MR207354L3V**

Serinc5 (NM_172588) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serinc5 (NM_172588) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Serinc5
Synonyms:	A130038L21Rik; AIGP3; C86123; TPO1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172588
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172588.2 , NP_766176.1
RefSeq Size:	5393 bp
RefSeq ORF:	1386 bp
Locus ID:	218442
UniProt ID:	Q8BHJ6
Cytogenetics:	13 C3



[View online »](#)

Gene Summary:

Restriction factor required to restrict infectivity of gammaretroviruses: acts by inhibiting early step of viral infection and impairing the ability of the viral particle to translocate its content to the cytoplasm (By similarity). Enhances the incorporation of serine into phosphatidylserine and sphingolipids. May play a role in providing serine molecules for the formation of myelin glycosphingolipids in oligodendrocytes (By similarity).[UniProtKB/Swiss-Prot Function]