

Product datasheet for **MR207360L4V**

Ints12 (NM_027927) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Ints12 (NM_027927) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Ints12 |
| Synonyms: | 1110020M19Rik; 2810027J24Rik; 4930529N21Rik; A230056J18Rik; A1666778; int12; Phf22 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_027927 |
| ORF Size: | 1386 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR207360). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_027927.3 |
| RefSeq Size: | 3014 bp |
| RefSeq ORF: | 1386 bp |
| Locus ID: | 71793 |
| UniProt ID: | Q9D168 |
| Cytogenetics: | 3 G3 |



[View online »](#)

Gene Summary:

Component of the Integrator complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes. Mediates recruitment of cytoplasmic dynein to the nuclear envelope, probably as component of the INT complex.[UniProtKB/Swiss-Prot Function]