

## Product datasheet for **MR207379L3V**

### Ctsc (NM\_009982) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctsc (NM_009982) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctsc
Synonyms:	AI047818; CatC; D; DP; DPP1; DPPI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009982
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207379).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009982.2</a>
RefSeq Size:	2472 bp
RefSeq ORF:	1389 bp
Locus ID:	13032
UniProt ID:	<a href="#">P97821</a>
Cytogenetics:	7 D3



[View online »](#)

**Gene Summary:**

This gene encodes a member of the peptidase C1 (papain) family of cysteine proteases. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate multiple protein products. These products include the dipeptidyl peptidase 1 light, heavy, and exclusion domain chains, which together comprise one subunit of the homotetrameric enzyme. This enzyme has amino dipeptidase activity and may play a role in the activation of granzymes during inflammation. Homozygous knockout mice for this gene exhibit impaired granzyme activation and enhanced survival in a sepsis model. [provided by RefSeq, Aug 2015]