

Product datasheet for **MR207385L4V**

Amotl2 (NM_019764) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Amotl2 (NM_019764) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Amotl2
Synonyms:	AW549739; C79691; Lccp; MASCOT; mKIAA0989
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019764
ORF Size:	2319 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207385).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019764.2 , NP_062738.2
RefSeq Size:	4191 bp
RefSeq ORF:	2319 bp
Locus ID:	56332
UniProt ID:	Q8K371
Cytogenetics:	9 F1



[View online »](#)

Gene Summary:

Regulates the translocation of phosphorylated SRC to peripheral cell-matrix adhesion sites. Required for proper architecture of actin filaments. Inhibits the Wnt/beta-catenin signaling pathway, probably by recruiting CTNNB1 to recycling endosomes and hence preventing its translocation to the nucleus. Participates in angiogenesis. May play a role in the polarity, proliferation and migration of endothelial cells. Selectively promotes FGF-induced MAPK activation through SRC (By similarity).[UniProtKB/Swiss-Prot Function]