

Product datasheet for **MR207403L3V**

2210010L05Rik (BC017534) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	2210010L05Rik (BC017534) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2210010L05Rik
Synonyms:	2210010L05Rik; 9630025I22Rik; AI836200; AW212394; MMR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC017534
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207403).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC017534 , AAH17534
RefSeq Size:	2744 bp
RefSeq ORF:	1391 bp
Locus ID:	98682
Cytogenetics:	1 C1.1
Gene Summary:	MHC class I receptor. Binds only to H-2 class I histocompatibility antigen, K-D alpha chain (H-2K(D)).[UniProtKB/Swiss-Prot Function]



[View online »](#)