

Product datasheet for **MR207422L4V**

Stk38 (NM_134115) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stk38 (NM_134115) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stk38
Synonyms:	5830476G13Rik; 9530097A09Rik; AA617404; Ndr1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134115
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207422).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134115.2
RefSeq Size:	3288 bp
RefSeq ORF:	1398 bp
Locus ID:	106504
UniProt ID:	Q91VJ4
Cytogenetics:	17 A3.3
Gene Summary:	Negative regulator of MAP3K1/2 signaling. Converts MAP3K2 from its phosphorylated form to its non-phosphorylated form and inhibits autophosphorylation of MAP3K2 (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)