

## Product datasheet for **MR207439L4V**

### Rnmt (NM\_026440) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnmt (NM_026440) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnmt
Synonyms:	2610002P10Rik; AI848273; mKIAA0398; Rg7mt1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026440
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207439).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026440.2</a> , <a href="#">NP_080716.1</a>
RefSeq Size:	5018 bp
RefSeq ORF:	1398 bp
Locus ID:	67897
UniProt ID:	<a href="#">Q9D0L8</a>
Cytogenetics:	18 E2
Gene Summary:	Catalytic subunit of the mRNA-capping methyltransferase RNMT:RAMAC complex that methylates the N7 position of the added guanosine to the 5'-cap structure of mRNAs. Binds RNA containing 5'-terminal GpppC.[UniProtKB/Swiss-Prot Function]



[View online »](#)