

Product datasheet for **MR207455L3V**

Scg3 (NM_001164790) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scg3 (NM_001164790) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scg3
Synonyms:	1B1075; AI385542; Chgd; SgIII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164790
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207455).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001164790.1 , NP_001158262.1
RefSeq Size:	2192 bp
RefSeq ORF:	1401 bp
Locus ID:	20255
UniProt ID:	P47867
Cytogenetics:	9 42.32 cM



[View online »](#)

Gene Summary:

Member of the granin protein family that regulates the biogenesis of secretory granules (PubMed:16219686). Acts as a sorting receptor for intragranular proteins including chromogranin A/CHGA (PubMed:12388744). May also play a role in angiogenesis (PubMed:28330905). Promotes endothelial proliferation, migration and tube formation through MEK/ERK signaling pathway (By similarity).[UniProtKB/Swiss-Prot Function]