

## Product datasheet for **MR207458L1V**

### Sox10 (NM\_011437) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sox10 (NM_011437) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sox10
Synonyms:	Dom; gt; Sox21
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_011437
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207458).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011437.1</a> , <a href="#">NP_035567.1</a>
RefSeq Size:	2713 bp
RefSeq ORF:	1401 bp
Locus ID:	20665
UniProt ID:	<a href="#">Q04888</a>
Cytogenetics:	15 37.7 cM



[View online »](#)

**Gene Summary:**

Transcription factor that plays a central role in developing and mature glia (PubMed:24204311, PubMed:27532821). Specifically activates expression of myelin genes, during oligodendrocyte (OL) maturation, such as DUSP15 and MYRF, thereby playing a central role in oligodendrocyte maturation and CNS myelination (PubMed:24204311, PubMed:27532821). Once induced, MYRF cooperates with SOX10 to implement the myelination program (PubMed:24204311). Transcriptional activator of MITF, acting synergistically with PAX3 (By similarity).[UniProtKB/Swiss-Prot Function]