

Product datasheet for **MR207473L4V**

Lmbrd1 (NM_026719) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lmbrd1 (NM_026719) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lmbrd1
Synonyms:	0910001K20Rik; AV347960
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026719
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207473).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026719.2
RefSeq Size:	5422 bp
RefSeq ORF:	1614 bp
Locus ID:	68421
UniProt ID:	Q8K0B2
Cytogenetics:	1 A5
Gene Summary:	Probable lysosomal cobalamin transporter. Required to export cobalamin from lysosomes allowing its conversion to cofactors (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)