

Product datasheet for **MR207488L4V**

Gtf2a1l (NM_023630) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtf2a1l (NM_023630) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gtf2a1l
Synonyms:	1700011N16Rik; Alf; Gtf2a1lf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023630
ORF Size:	1407 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207488).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023630.1 , NP_076119.2
RefSeq Size:	1619 bp
RefSeq ORF:	1407 bp
Locus ID:	71828
UniProt ID:	Q8R4I4
Cytogenetics:	17 E4
Gene Summary:	May function as a testis specific transcription factor. Binds DNA in conjunction with GTF2A2 and TBP (the TATA-binding protein) and together with GTF2A2, allows mRNA transcription. [UniProtKB/Swiss-Prot Function]



[View online »](#)