

Product datasheet for MR207512L4

Rufy3 (NM_027530) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rufy3 (NM_027530) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Rufy3
Synonyms:	2810428M05Rik; 6330416M07Rik; AW455998; AW538594; D5Bwg0860e; Ripx; Rpipx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207512).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_027530
ORF Size:	1410 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027530.2 , NP_081806.1
RefSeq Size:	4097 bp
RefSeq ORF:	1410 bp
Locus ID:	52822
UniProt ID:	Q9D394
Cytogenetics:	5 43.77 cM
Gene Summary:	Plays a role in the generation of neuronal polarity formation and axon growth (PubMed:24720729). Implicated in the formation of a single axon by developing neurons (PubMed:24720729). May inhibit the formation of additional axons by inhibition of PI3K in minor neuronal processes (By similarity). Plays a role in the formation of F-actin-enriched protrusive structures at the cell periphery (By similarity). Plays a role in cytoskeletal organization by regulating the subcellular localization of FSCN1 and DBN1 at axonal growth cones (PubMed:24720729). Promotes gastric cancer cell migration and invasion in a PAK1-dependent manner (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR207512L4