

Product datasheet for **MR207515L3V**

Snx17 (NM_153680) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snx17 (NM_153680) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snx17
Synonyms:	5830447M19Rik; A1790646; b2b1625.1Clo; D5Ert260e; mK1AA0064
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153680
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207515).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153680.2
RefSeq Size:	1941 bp
RefSeq ORF:	1413 bp
Locus ID:	266781
UniProt ID:	Q8BVL3
Cytogenetics:	5 B1



[View online »](#)

Gene Summary:

Critical regulator of endosomal recycling of numerous surface proteins, including integrins, signaling receptor and channels (PubMed:12169628, PubMed:16052210). Binds to NPxY sequences in the cytoplasmic tails of target cargos (By similarity). Associates with retriever and CCC complexes to prevent lysosomal degradation and promote cell surface recycling of numerous cargos such as integrins ITGB1, ITGB5 and their associated alpha subunits (By similarity). Also required for maintenance of normal cell surface levels of APP and LRP1 (PubMed:16052210, PubMed:18276590). Interacts with membranes containing phosphatidylinositol 3-phosphate (PtdIns(3P)) (By similarity).[UniProtKB/Swiss-Prot Function]