

## Product datasheet for **MR207558L4V**

### **Ptdss2 (NM\_013782) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ptdss2 (NM_013782) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ptdss2
Synonyms:	AI481134; AI988017; PSS2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013782
ORF Size:	1419 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207558).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013782.4</a> , <a href="#">NP_038810.2</a>
RefSeq Size:	2172 bp
RefSeq ORF:	1422 bp
Locus ID:	27388
UniProt ID:	<a href="#">Q9Z1X2</a>
Cytogenetics:	7 F5



[View online »](#)

**Gene Summary:**

Catalyzes a base-exchange reaction in which the polar head group of phosphatidylethanolamine (PE) or phosphatidylcholine (PC) is replaced by L-serine. PTDSS2 is specific for phosphatidylethanolamine and does not act on phosphatidylcholine.  
[UniProtKB/Swiss-Prot Function]