

Product datasheet for **MR207578L2V**

Sptlc1 (NM_009269) Mouse Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	Sptlc1 (NM_009269) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sptlc1
Synonyms:	AW552086; C77762; E030036H05; Lcb1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_009269
ORF Size:	1419 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207578).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009269.2 , NP_033295.2
RefSeq Size:	2601 bp
RefSeq ORF:	1422 bp
Locus ID:	268656
UniProt ID:	O35704
Cytogenetics:	13 B1



[View online »](#)

Gene Summary:

Serine palmitoyltransferase (SPT) (PubMed:28100772). The heterodimer formed with SPTLC2 or SPTLC3 constitutes the catalytic core. The composition of the serine palmitoyltransferase (SPT) complex determines the substrate preference. The SPTLC1-SPTLC2-SPTSSA complex shows a strong preference for C16-CoA substrate, while the SPTLC1-SPTLC3-SPTSSA isozyme uses both C14-CoA and C16-CoA as substrates. The SPTLC1-SPTLC2-SPTSSB complex displays a strong preference for C18-CoA substrate, while the SPTLC1-SPTLC3-SPTSSB isozyme has the ability to use a broader range of acyl-CoAs (By similarity). Required for adipocyte cell viability and metabolic homeostasis (PubMed:28100772).[UniProtKB/Swiss-Prot Function]

Product images:

[MR207578L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR207578L2V particle to overexpress human Sptlc1-mGFP fusion protein.