

## Product datasheet for **MR207597L4V**

### **Bpifb1 (NM\_001012392) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bpifb1 (NM_001012392) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bpifb1
Synonyms:	Lplunc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001012392
ORF Size:	1425 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207597).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001012392.1</a> , <a href="#">NP_001012392.1</a>
RefSeq Size:	1690 bp
RefSeq ORF:	1425 bp
Locus ID:	228801
UniProt ID:	<a href="#">Q61114</a>
Cytogenetics:	2 H1



[View online »](#)

**Gene Summary:**

May play a role in innate immunity in mouth, nose and lungs. Binds bacterial lipopolysaccharide (LPS) and modulates the cellular responses to LPS. May be involved in formation of the left-right axis in the node of the developing embryo.[UniProtKB/Swiss-Prot Function]