

## Product datasheet for **MR207634L4V**

### Gc (NM\_008096) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gc (NM_008096) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gc
Synonyms:	DBP; VDB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008096
ORF Size:	1428 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207634).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008096.1</a>
RefSeq Size:	1681 bp
RefSeq ORF:	1431 bp
Locus ID:	14473
UniProt ID:	<a href="#">P21614</a>
Cytogenetics:	5 44.32 cM
Gene Summary:	Involved in vitamin D transport and storage, scavenging of extracellular G-actin, enhancement of the chemotactic activity of C5 alpha for neutrophils in inflammation and macrophage activation.[UniProtKB/Swiss-Prot Function]



[View online »](#)