

Product datasheet for **MR207652L3V**

Micu1 (NM_144822) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Micu1 (NM_144822) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Micu1
Synonyms:	C730016L05Rik; Calc; Cbara1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144822
ORF Size:	1434 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207652).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144822.2
RefSeq Size:	2358 bp
RefSeq ORF:	1434 bp
Locus ID:	216001
UniProt ID:	Q8VCX5
Cytogenetics:	10 B4



[View online »](#)

Gene Summary:

Key regulator of mitochondrial calcium uniporter (MCU) that senses calcium level via its EF-hand domains (PubMed:24560927). MICU1 and MICU2 form a disulfide-linked heterodimer that stimulates and inhibits MCU activity, depending on the concentration of calcium (PubMed:24560927). MICU1 acts both as an activator or inhibitor of mitochondrial calcium uptake (By similarity). Acts as a gatekeeper of MCU at low concentration of calcium, preventing channel opening (By similarity). Enhances MCU opening at high calcium concentration, allowing a rapid response of mitochondria to calcium signals generated in the cytoplasm (PubMed:24560927). Regulates glucose-dependent insulin secretion in pancreatic beta-cells by regulating mitochondrial calcium uptake (By similarity). Induces T-helper 1-mediated autoreactivity, which is accompanied by the release of IFNG (By similarity). [UniProtKB/Swiss-Prot Function]