

Product datasheet for **MR207655L3V**

Ppox (NM_008911) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppox (NM_008911) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppox
Synonyms:	Ppo
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008911
ORF Size:	1434 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207655).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008911.2 , NP_032937.1
RefSeq Size:	1821 bp
RefSeq ORF:	1434 bp
Locus ID:	19044
UniProt ID:	P51175
Cytogenetics:	1 H3
Gene Summary:	Catalyzes the 6-electron oxidation of protoporphyrinogen-IX to form protoporphyrin-IX. [UniProtKB/Swiss-Prot Function]



[View online »](#)