

## Product datasheet for **MR207663L4V**

### Shmt1 (NM\_009171) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Shmt1 (NM_009171) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Shmt1
Synonyms:	AI324848; AI385541; C81125; mshmt; mshmt1; mshmt2; Shmt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009171
ORF Size:	1434 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207663).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009171.2</a> , <a href="#">NP_033197.2</a>
RefSeq Size:	2074 bp
RefSeq ORF:	1437 bp
Locus ID:	20425
UniProt ID:	<a href="#">P50431</a>
Cytogenetics:	11 B2
Gene Summary:	Interconversion of serine and glycine.[UniProtKB/Swiss-Prot Function]



[View online »](#)