

Product datasheet for **MR207686L4V**

Irak2 (BC085324) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Irak2 (BC085324) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Irak2
Synonyms:	6330415L08Rik; AI649099; IRAK-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC085324
ORF Size:	1437 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207686).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC085324 , AAH85324
RefSeq Size:	3049 bp
RefSeq ORF:	1439 bp
Locus ID:	108960
Cytogenetics:	6 52.82 cM
Gene Summary:	Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.[UniProtKB/Swiss-Prot Function]



[View online »](#)